

The logo for Delta⁴ by ScandiDos. The word "Delta" is in a bold, dark red font, and the superscript "4" is in a lighter red. Below it, "by ScandiDos" is written in a smaller, grey font.

Delta⁴
by ScandiDos

Delta⁴ TPV

A unique Monte Carlo based Treatment
Planning Verification in 4D

What is Delta⁴ TPV?

- Secondary dose calculation tool
- With a Monte Carlo dose engine
- Fast *and* accurate

Why Monte Carlo?

- Gold standard
 - *Monte Carlo is the gold standard computation method for accurate dose calculations in radiotherapy*

Why Monte Carlo?

- Speed
 - *Monte Carlo computation time scales very well with the number of cores*

Why Monte Carlo?

- ICRU 83 recommendations
 - *“Therefore, as an alternative to a set of measured absorbed-dose distributions, it is acceptable to use an **independent absorbed-dose calculation that is at least as accurate** as the absorbed-dose calculation being tested as previously verified against commissioning measurements.”*

Cooperation with Radyalis

Radyalis

- Exclusive cooperation with American company Radyalis
- Radyalis are experts on Monte Carlo calculations for radiotherapy applications

What is Delta⁴ TPV?

Dose QA

and

Plan QA

- Target covered?
- OAR dose acceptable?
 - RTOG protocols

Where does it fit in the Delta⁴ Family?



Plan Verification

- Delta⁴ TPV

Pre-treatment Verification

- Delta⁴ Phantom+

At-treatment Verification

- Delta⁴ Discover

Good plan...

delivered correctly...

every day.

What is the Delta⁴ TPV workflow?

Automated Process

- Export DICOM files from TPS
- Delta⁴ TPV automatically
 - Starts the Monte Carlo calculation
 - Evaluates the results
 - Saves the results in the Delta⁴ database

What is the Delta⁴ TPV workflow?

Evaluation

- Overview




What is the Delta⁴ TPV workflow?

Deeper evaluation

- Gamma
- DVH
- Plan objectives




Structure	Objective	TPS	TPV
PTV45 and PTV50	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.10$	✘ 0.09	✔ 0.10
PTV50_krans10mm	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.10$	✔ 0.35	✔ 0.36
PTV50_krans5mm	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.10$	✔ 0.20	✔ 0.23
Bladder	$D_{70\%} \leq 40.0 \text{ Gy}$	✘ 41.2 Gy	✘ 41.2 Gy
	$D_{55\%} \leq 50.0 \text{ Gy}$	✔ 47.4 Gy	✔ 47.4 Gy
	$D_{30\%} \leq 70.0 \text{ Gy}$	✔ 50.0 Gy	✔ 50.0 Gy
	$D_{15\%} \leq 80.0 \text{ Gy}$	✔ 50.5 Gy	✔ 50.7 Gy
Kidney L	$D_{100\%} \leq 23.0 \text{ Gy}$	✔ 2.9 Gy	✔ 2.9 Gy
	$D_{67\%} \leq 30.0 \text{ Gy}$	✔ 5.7 Gy	✔ 5.9 Gy
	$D_{33\%} \leq 50.0 \text{ Gy}$	✔ 9.2 Gy	✔ 9.3 Gy




What is the Delta⁴ TPV workflow?

Approval



What is the Delta⁴ TPV workflow?

Target	OAR	Gamma	Approval
Structure	Objective	TPS	TPV
GTV-Ln	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.01$	0,04	0,04
GTV-T	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.01$	0,05	0,05
PTV45 and PTV50	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.01$	0,09	0,10
PTV50_krans10mm	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.01$	0,35	0,37
PTV50_krans5mm	$HI = (D_{2\%} - D_{98\%}) / D_{50\%} \geq 0.01$	0,20	0,23
Bladder	$D_{70\%} \leq 40.0 \text{ Gy}$	✗ 41,2 Gy	✗ 41,0 Gy
	$D_{55\%} \leq 50.0 \text{ Gy}$	✓ 47,4 Gy	✓ 46,9 Gy
	$D_{30\%} \leq 70.0 \text{ Gy}$	✓ 50,0 Gy	✓ 49,9 Gy
	$D_{15\%} \leq 80.0 \text{ Gy}$	✓ 50,5 Gy	✓ 50,6 Gy
Kidney L	$D_{100\%} \leq 23.0 \text{ Gy}$	✓ 2,9 Gy	✓ 2,9 Gy
	$D_{67\%} \leq 30.0 \text{ Gy}$	✓ 5,7 Gy	✓ 5,9 Gy
	$D_{33\%} \leq 50.0 \text{ Gy}$	✓ 9,2 Gy	✓ 9,3 Gy
Kidney R	$D_{100\%} \leq 23.0 \text{ Gy}$	✓ 5,9 Gy	✓ 5,6 Gy
	$D_{67\%} \leq 30.0 \text{ Gy}$	✓ 8,5 Gy	✓ 8,8 Gy
	$D_{33\%} \leq 50.0 \text{ Gy}$	✓ 10,7 Gy	✓ 10,8 Gy
Rektum	$D_{200 \text{ cc}} \leq 30.0 \text{ Gy}$?	?
	$D_{150 \text{ cc}} \leq 35.0 \text{ Gy}$?	?
	$D_{100 \text{ cc}} \leq 40.0 \text{ Gy}$?	?
	$D_{20 \text{ cc}} \leq 45.0 \text{ Gy}$	✗ 48,8 Gy	✗ 47,9 Gy
	$D_{\text{Max}} \leq 52.0 \text{ Gy}$	✓ 51,9 Gy	✓ 51,3 Gy
Sigmoid	$D_{\text{Max}} \leq 5.0 \text{ Gy}$	✗ 51,6 Gy	✗ 52,0 Gy
Spinal cord	$D_{\text{Max}} \leq 45.0 \text{ Gy}$	✓ 38,4 Gy	✓ 37,8 Gy



What hardware do I need?

TPV Calculation server

–Powerful PC

Delta⁴ Viewing station

–Evaluation and Approval

–Regular PC

What hardware do I need?

TPV Calculation server

- Runs on any PC with Windows 10
- For speed use at least
 - Processor: Intel® Xeon® W with 18 cores/36 threads/2.3GHz
 - RAM: 64 GB

How fast is Delta⁴ Monte Carlo?

*100 - 1000 times faster than general
Monte Carlo codes*

How fast is Delta⁴ Monte Carlo?

Plan Type	Body Site	Energy Machine Type	Time (seconds)
VMAT	H&N	6 MV Elekta Agility	371
VMAT	Thorax	6 MV Elekta Agility	592
VMAT	Cervix	6 MV Varian TrueBeam Millennium	715
3D conformal	Spine	10 MV Varian TrueBeam Millennium	102
VMAT	Prostate	6 MV FFF Varian TrueBeam Millennium	318

How accurate is Delta⁴ Monte Carlo?

Clinical test cases

- Measurements
- Gamma 2%/2mm: Pass-rate 96%
For voxels with more than 10% of dose max

Physics test cases

- EGSnrc
- Geant4
- Gamma 1%/1mm: Pass-rate 96%
For voxels with more than 10% of dose max

Why is Delta⁴ Monte Carlo fast *and* accurate?

- Optimized for radiotherapy applications
- Efficient algorithms, modeling, and data representation
- Utilizes all available computing resources

Why is Delta⁴ Monte Carlo fast *and* accurate?

- Simulates all relevant particles and materials
- Includes all relevant physics, models arbitrary geometries
- Proprietary methods for convergence, statistical noise controls and dose uncertainty bounds

How is commissioning done?

- Same measurements as for TPS commissioning
- Commissioning done by ScandiDos
- Graphical tools for evaluation of the result

Delta⁴

by ScandiDos

Secondary dose calculation tool
- within the Delta⁴ family or Standalone

Time efficient with Automated Process
Highest Accuracy with Monte Carlo
Effective and Objective Evaluation (RTOG)

Delta⁴ TPV

A unique Monte Carlo based Treatment
Planning Verification in 4D