

Maximize the performance of your Delta4 system. Keep your service contract active and your software up to date to benefit from the latest features and improvements, with additional capabilities available through optional licensed modules.

- Compensate for temperature variations in the detectors*



**The Automatic Temperature Correction function requires an up-to-date software and hardware service contract or a purchased upgrade. Not available for all systems. Please contact your ScandiDos representative for further information.

See How Delta4 Software Has Evolved

Delta4 software is continuously developed with new functionality and workflow improvements.

See what has been added since your current software version.

NOV 2021

Daily Position Corrections

Automatically apply position offsets for plans measured on the same day.

Test Data Management

Save plans as new patients for easier management of calibration and test data.

Synchrony kV Dose Separation

Prevent imaging dose from distorting measurements by automatically separating kV dose.

MAR 2023

Merge Multiple Measurements

No limitation in longitudinal field size and resolution for measurements.

Machine QA Support for Halcyon

MAY 2024

Statistics Module

View trends, track changes and identify outliers in measurement data.

Automatic Temperature Correction

Automatically measure the temperature in each individual diode at the start of each measurement.

JUN 2026

Machine QA 3D

Use a highly modulated clinical plan to monitor overall linac performance over time.

Machine QA Output Check

Monitor linearity, output constancy, and output at different dose rates.

Statistics Module Improvements

Use linac tagging to track and compare individual linacs.

LATEST VERSION

OCT 2024

Statistics Module Improvements

Mine the database, create comparison diagrams, and copy plans for comparisons.

MAY 2025

Machine QA Protocols

Optimize workflow with tailored QA protocols. Deliver all checks in a single DICOM plan.

Improved Merge Function

Merge measurements in normal and reversed phantom positions to avoid irradiation of electronics.

NOV 2025

Machine QA Leaf Leakage Check

Monitor intraleaf leakage over time.



Get the Latest from Your Delta4 System

Update to the latest software version to access the newest features and improvements available for your system and licensed modules.

✉ support@scandidos.com

☎ +46 18 472 30 30

Keep Your Delta4 System Performing at Its Best

Get more from your Delta4 system with comprehensive service and support. A ScandiDos service contract provides access to software and firmware updates, training, technical expertise, and priority support.

With a service contract in place, you can keep your Delta4 system current as new software and firmware updates become available.



EXPERIENCE PEACE OF MIND - ORDER YOUR SERVICE CONTRACT TODAY

Corporate Headquarters

ScandiDos AB (publ)
Uppsala Science Park
SE-751 83 Uppsala, Sweden
+46 18 472 30 30
info@scandidos.com
delta4family.com
 scandidos

© 2026 ScandiDos AB

To Order Your Service Contract

Global Support Team
support@scandidos.com
+46 18 472 30 30

USA Service and Support
+1 608 630 8393

France Service and Support
+33 6 31 50 45 81